DK50

DK50-IO DK50 PLUS

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IMPORTANT INFORMATION

1. CE MARKING

Products labeled with the CE mark of compliance meet the safety guidelines (93/42/EEC) of the European Union.

2. WARNINGS

2.1. General warnings

- This Installation, Operation and Maintenance Manual is a part of the appliance and must be kept with the compressor. Careful review of this manual will provide the information necessary for correct operation of the appliance.
- The safety of operating personnel and trouble-free operation of the appliance are guaranteed only if original parts are used. Only accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other accessories or consumable materials are used, the manufacturer cannot be held responsible for the safe operation of the appliance. This guarantee does not cover damages originating from the use of accessories or consumable material other than those specified or suggested by the manufacturer.
- The manufacturer guarantees the safety, reliability and function of the appliance only if:
 - Installation, new settings, amendments, extensions and repairs are performed by the manufacturer or its representative, or a service provider authorized by the manufacturer
 - The appliance is used in accordance with this Installation, Operation and Maintenance Manual
- The manufacturer reserves all rights for the protection of its wiring diagrams, methods and names.
- Translation of Manual for Installation, Operation and Maintenance is carried out in accordance with the best knowledge. In the case of ambiguities, the Slovak version of the text prevails.

2.2. General safety warnings

The manufacturer developed and designed the equipment in such a way so that any risks were excluded if it is used according to intention. The manufacturer considers it to be its obligation to describe the following safety measures in order to exclude residual damages.

- Operation of the appliance must be in compliance with all local codes and regulations.
- Original packaging should be kept for the return of the appliance. Only the original packaging ensures protection of the appliance during transport. If it is necessary to return the appliance during the guarantee period, the manufacturer is not liable for damages caused by improper packaging.
- Each time the appliance is used, the operator must make sure that it is functioning correctly and safely.
- The user must fully understand the operation of the appliance.
- The product is not intended for operation in areas with a risk of explosion.
- If any problem occurs during use of the appliance, the user must inform his supplier immediately.

2.3. Electrical system safety warnings

- The appliance must be connected to earth (grounded).
- Before the appliance is plugged in, make sure that the mains voltage and mains frequency stated on the appliance are the same as the power mains.
- Prior to putting into operation it is necessary to check for possible damage of the equipment and connected air and electric distributions. Damaged pneumatic and electric lines must be immediately replaced.
- Immediately disconnect the appliance from the mains (pull out mains plug) if a technical failure occurs.
- During repairs and maintenance, ensure that:
 - The mains plug is pulled out from the socket
 - Pressure pipes are vented and pressure is released from the air tank.
- The appliance must be installed by an approved, qualified technician.

3. ALERT NOTICES AND SYMBOLS

In the Installation, Operation and Maintenance Manual and on the appliance and its packaging, the following labels or symbols are used for important information:

	Information, instructions and cautions for the prevention of damage to health or materials
	Caution! Dangerous electric voltage
	Information regarding correct use of the appliance and other warnings
	CE mark of compliance
	Caution! Hot surface
	Compressor is remote-controlled and may start without warning
	Earth (ground) connection
	Terminal for ground connection
	Fuse
	Alternating current
	Handling mark on package – Fragile, handle with care
	Handling mark on package – This way up (vertical position of cargo)
	Handling mark on package – Protect against moisture
	Handling mark on package – Temperature during storage and transport
	Handling mark on package – Limited stacking
	Mark on package – Recyclable material

4. STORAGE AND TRANSPORT

The compressor is shipped in cardboard that protects the appliance from damage during transport.

Caution! For transport, always use the original packaging and secure the compressor in the upright position.

Protect the compressor from humidity and extreme temperatures during transport and storage. A compressor in its original packaging can be stored in a warm, dry and dust-free area. Do not store near any chemical substances.

Keep packaging material if possible. If not, please dispose of the packaging material in an environmentally friendly way and recycle if possible.

Caution! Before moving or transporting the compressor, release all the air pressure from the tank and hoses and drain the condensed water.

5. TECHNICAL DATA

	DK50 Z	DK50 S	DK50-10 Z	DK50-10 S	DK50 PLUS	DK50 PLUS S
Nominal voltage / frequency (*) V / Hz	230 / 50 230 / 60 110 / 60	230 / 50 230 / 60 110 / 60	230 / 50 230 / 60 110 / 60	230 / 50 230 / 60 110 / 60	230 / 50 230 / 60 110 / 60	230 / 50 230 / 60 110 / 60
Efficiency of compressor at over-pressure 5 bar Lit.min ⁻¹	75	75	75	75	75	75
Efficiency of compressor with dryer at over-pressure 5 bar Lit.min ⁻¹	-	-	75	75	75	75
Efficiency of compressor with KJF-1 at over-pressure 5 bar Lit.min ⁻¹	75	75	75	75	75	75
Maximal current A	3.4 4.3 8.6	3.4 4.3 8.6	3.4 4.3 8.6	3.4 4.3 8.6	3.4 4.3 8.6	3.4 4.3 8.6
Maximal current of compressor with dryer A	-	-	3.6 4.5 8.8	3.6 4.5 8.8	3.6 4.5 8.8	3.6 4.5 8.8
Motor performance kW	0.55	0.55	0.55	0.55	0.55	0.55
Air tank capacity Lit.	5	5	10	10	25	25
Pressure range (**) bar	4.5 – 6.0	4.5 – 6.0	4.5 – 6.0	4.5 – 6.0	4.5 – 6.0	4.5 – 6.0
Maximum operating pressure of safety valve bar	8.0	8.0	8.0	8.0	8.0	8.0
Sound level L _{ptA} [dB]	≤ 65	≤ 45	≤ 65	≤ 45	≤ 66	≤ 47
Mode of operation	continual S 1	continual S 1	continual S 1	continual S 1	continual S 1	intermittent S 3 – 80%
Mode of operation of compressor with dryer	-	-	intermittent S 3 – 50 %	intermittent S 3 – 50%	intermittent S 3 – 60%	intermittent S 3 – 60%
Dimensions of compressor / of compressor with dryer W x L x H mm	290x430x490 /-	380x525x575 /-	330x430x530 / 330x580x570	420x525x620 / 420x675x620	460x460x690 / 460x515x690	560x665x860
Weight of compressor / of compressor with dryer kg	34/ -	46/-	36/47	49/64	48 / 58	83 / 93
Weight of compressor / of compressor with dryer Packaget in cardboard kg	36/ -	48/ -	38/50	52/67	52 / 62	90 / 98
Drying point of compressor Atmospheric condensation point	-	-	to -20°C	to -20°C	to -20°C	to -20°C
Version EN 60 601-1	Type B, class I.					

Notices:

- * When ordering, state the version of compressor
- ** Range of pressure: consult with contractor
- Weight of compressor with KJF-1 increase about 3kg

Climatic conditions during storage and transport

Temperature -25°C to +55°C, 24 h to +70°C
Relative air humidity 10% to 90% (no condensation)

Climatic operation conditions

Temperature +5°C to +40°C
Relative air humidity 70%

6. PRODUCT DESCRIPTION

6.1. Model variations and their uses

Compressors are the source of clean, oil-free compressed air used to drive dental appliances and equipment.

Compressors models are designed for the following uses:

Dental compressors DK50 Z and DK50-10 Z - sit on a free-standing base.

Dental compressors DK50 Z/K and DK 50-10 Z/K - sit on a free-standing base and feature a condensation and filtration unit (KJF1.)

Dental compressors DK50-10 Z/M - sit on a free-standing base and feature an adsorption air dryer M1.10.

Dental compressors DK50 S and DK50-10 S - feature compact soundproof boxes suitable for placing in a dentist's office.

Dental compressors DK50 S/K and DK50-10 S/K - feature compact boxes and a condensation and filtration unit (KJF1).

Dental compressors DK50-10 S/M - feature compact boxes and an adsorption air dryer M1.10.

Dental compressors DK50 PLUS - sit on a free-standing base

Dental compressors DK50 PLUS/K - sit on a free-standing base and feature a condensation and filtration unit (KJF1.)

Dental compressors DK50 PLUS/M - sit on a free-standing base and feature an adsorption air dryer M1a.

Dental compressors DK50 PLUS S - feature compact soundproof boxes suitable for placing in a dentist's office.

Dental compressors DK50 PLUS S/K - feature compact boxes and a condensation and filtration unit (KJF1).

Dental compressors DK50 PLUS S/M - feature compact boxes and an adsorption air dryer M1a.

DK50 Z

DK50-10 Z

DK50-10 S

DK50 PLUS

DK50 PLUS S

KJF1

MONZUN
M1.10MONZUN
M1a

Without additional filtration equipment, the compressed air from a compressor is not suitable for the operation of breathing appliances or similar equipment.

7. FUNCTION

Compressor (Fig.1)

The compressor (1) draws in air through a filter (8) and compresses it through a check valve (3) into an air tank (2). The connected apparatus draws the compressed air from the air tank until the pressure drops to a default preset level on the air-pressure switch (4) switching the compressor on. The compressor again compresses air into the nozzle until the maximum pressure is reached and the compressor switches off. After compressor aggregate is switched off, pressure hose shall be pressure-release solenoid valve (13). Safety valve (5) prevents the pressure in air chamber from rising above the maximal allowed value. The drain valve (7) releases the condensate from the air nozzle. Compressed, clean air free from oil traces is stored in the air tank ready for use.

Compressor with dryer M1.10 (Fig.2, Fig.4)

The compressor (1) draws in air through a filter (8) and sends it to the air dryer in compressed form. The air proceeds through the cooler (15) and through the dryer chamber (9) with dessicant (16). The moisture is captured in a output filter (17) and check valve (3), sending clean, dry air into the air tank (2). Adsorber shall be regenerated after every switching off of compressor using pressure switch during dryer chamber depressurisation. Air escapes from adsorption chamber via open solenoid valve (14), while it is concurrently blown with dried air. The dessicant is regenerated and the captured water from the chamber flows out through the drain valve. Dry, clean compressed air free from oil traces is stored in the air tank ready for use.

Additional dryer regeneration - without case (Valid for M1.10)

To ensure efficient operation of the dryer, the dessicant must be regenerated twice a month. Open the valve for additional regeneration (19) for 8 hours. During this time the dessicant shall regenerate, ensuring its high efficiency for the life of the dryer. After this time of forced regeneration it is necessary to close the regeneration valve by (fine) clockwise rotation.

Additional dryer regeneration - in a case (Valid for M1.10) (Fig.5)

For the compressor with dryer MONZUN M1.10 in a case, the valve operation for regeneration (19) connects through a flexible spindle to the front panel of the case. To ensure efficient operation of the dryer, the dessicant must be regenerated twice a month. Open the valve for additional regeneration by turning the button (21) from position "0" to position "I" to the left (counterclockwise) for 8 hours. After the regeneration is complete, close the regeneration valve by turning the button from position "I" to position "0" (clockwise).

Compressor with dryer M1a (Fig.3)

The compressor (1) draws in air through a filter (8) and sends it to the air dryer in compressed form. The air proceeds through the cooler (15) and through the dryer chamber (9) with dessicant (16). The moisture is captured in a output filter (17) and check valve (3), sending clean, dry air into the air tank (2). Adsorber shall be regenerated after every switching off of compressor using pressure switch during dryer chamber depressurisation. Air escapes from adsorption chamber via open solenoid valve (14), while it is concurrently blown with dried air. The level of drying of the passing-through air is monitored by a humidity sensor (23), that shall open the regeneration valve (24) in the case of registering humidity exceeding the set humidity and at the time of a pause of the compressor it shall automatically regenerate the drying agent in the chamber. The entrapped water from the chamber shall be pushed out via a valve. Dry, clean compressed air free from oil traces is stored in the air tank ready for use.

Compressor with condensation and filtration unit (Fig.7)

The compressor (1) draws in air through a filter (8) and compresses it through a check valve (3) into an air tank (2). The compressed air from the nozzle flows through a cooler (10) that cools the compressed air. The condensed moisture is trapped in the filter (11) and automatically separates as condensate (12). Dried, clean compressed air, free from oil traces, is ready for use.

Compressor box (Fig.5, Fig.6)

The soundproof box is compact yet allows sufficient exchange of cooling air. It can be placed in a dentist's office. The ventilator under the aggregate of a compressor provides cooling of compressor and it is in operation at the same time with an engine of the compressor. After prolonged use the temperature in the case may rise above 40°C, causing the cooling fan blower to automatically turn on. After cooling the case area to 32°C the fan blower turns off automatically. Door of the casing with right opening may be changed to left opening (see Chapter 9).


Make sure that nothing impedes the free flow of air under and around the compressor. Never cover the hot air outlet on the top back side of the case.

If placing the compressor on a soft floor such as carpet, create space for ventilation between the base and floor or the box and floor, e.g. underpin the footings with hard pads.

Fig.1 – Compressor

1. Compressor motor
2. Air tank
3. Check valve
4. Pressure switch
5. Safety valve
6. Manometer
7. Drain valve
8. Input filter
9. Dryer chamber
10. Pipe cooler
11. Output filter
12. Condenser outlet
13. Solenoid valve
14. Solenoid valve of dryer
15. Cooler
16. Dessicant
17. Output filter
18. Sieve
19. Valve for regeneration
20. Stopper
21. Regeneration valve button
22. Bottle
23. Humidity sensor
24. Regeneration valve
25. Box
26. Lock
27. Connecting reinforcement
28. Wall stopper
29. Switch
30. Manometer
31. Magnetic bottle holder
32. Door hinge
33. Wheels
34. Socket on the box
35. Silencer
40. Hose of manometer
41. Compressor Fan
43. Output valve
44. Rectification screw
45. Door pin
46. Compressor handle

Fig. 2 - Compressor with air dryer MONZUN-M1.10


Fig. 3 - Compressor with air dryer MONZUN- M1a


Fig.4 - Air dryer MONZUN-M1.10

Fig.5 – Compressor in box


Fig.6 - Compressor in box DK50 PLUS S



Fig.7 - Compressor with condensation and filtration unit KJF1

INSTALLATION

8. USE

- Aspirators (suction unit) are designed for operation in dry, well ventilated and dust-free area where ambient temperature is within the range of +5°C to +40°C and relative air humidity does not exceed 70%. The compressor must be installed so that it is accessible at all times for operating and maintenance. Please ensure that the appliance label is accessible as well.
- The appliance must stand on a flat, sufficiently stable base. See paragraph 5 (Technical data) when positioning or lifting the compressor.
- Compressors cannot be exposed to outdoor environments. The appliance cannot be used in moist or wet environments. Do not use the compressor in the presence of explosive gases, dust or combustible liquids.
- Before connecting the compressor to medical equipment, the supplier must confirm that it meets all requirements for its use. Refer to the technical data of the product for this purpose. When a unit is to be built-in, classification and evaluation of compatibility must be done by the manufacturer or supplier of the product to be used.
- Any use other than that described in this manual is not covered by the guarantee, and the manufacturer is not liable for any damages that may result. The operator/user assumes all risk.

9. INSTALLATION

Only qualified personnel can install and start up the appliance and train operating personnel in its correct use and maintenance. Installation and training of all operators shall be confirmed by the installer's signature on the certificate of installation.


Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.

Caution! When in operation, the compressor is hot. Burns or fire may result if contact is made by the operator or any flammable material.

9.1. Placement of the compressor

Dental compressor with base DK50 Z, DK50-10 Z, DK50-10 Z/M, DK50 PLUS, DK50 PLUS/M (Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Direct the output pressure hose, drain hose and power cord out the back of the compressor.

Dental compressor in box DK50 S, DK50-10 S (Fig.5, Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Direct the output pressure hose, drain hose and power cord out the back of the compressor. Slide the box over the compressor so that the front face of the box matches the front part of the compressor and the box is fully seated. Make sure that the pressure hose, drain hose and electric cord come out via the opening at the back of the box. Position the drain hose with its valve in the holder at the rear of the box.

Dental compressor in box DK50-10 S/M (Fig.5, Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Direct the output pressure hose, drain hose and power cord out the back of the compressor. Fit the housing over the top of the compressor, connect the flexible shaft to the control button, fasten with the screw and put the lid on the cabinet housing (see picture). Make sure that the pressure hose, drain hose and electric cord come out via the opening at the back of the box. Position the drain hose with its valve in the holder at the rear of the box.

Dental compressor in box DK50 PLUS S (Fig.6, Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the wall-mounted stopper (28) - 2 pcs onto the compressor housing in the rear top part of the housing and put the housing onto a required place. The stoppers provide a sufficient distance of the housing from a wall for thorough ventilation. For setting up the compressor in housing you must open the door on the housing using the attached key and remove connecting reinforcement (27) in the front bottom part of housing. If necessary, the door may be disassembled using door hinge (32). Protrude pressure hose via a hole in housing and connect it to the appliance in a suitable way. Grasp the compressor at its handle and put it into the casing using built-in wheels (33). Fix the hose of a manometer in a hose into the fast-on coupling on a compressor, put the connecting reinforcement (27) back and connect the pressure hose to a compressor. Insert the electric power cord of a compressor into a socket (34) on a housing. By slight rotation of rectification screws (44) set the correct position of door against the casing frame. When closing the door the pin (45) on the door must easily snap in the opening in the casing frame. Close the housing doors and duly lock the lock (26). Connect the mains plug into the mains socket.

It is not allowed to leave the key in a lock! It must be saved against non-instructed persons!

Dental compressor in box DK50 PLUS S/M (Fig.6, Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the compressor into housing similarly as in the previous paragraph. Prior to placing the compressor into housing, protrude house for condensate drain via hole in housing and connect it to a bottle (22). The magnetic holder (31) with a vessel (22), for entrapping condensate from a dryer may be fixed at the sides of housing or from the front on its doors. When fixing the holder with a vessel at the housing side it is necessary to consider a space of at least 11 cm between the housing and furniture. Distance smaller than the specified one may cause problem with handling of the vessel.

9.2. Compressed air outlet

(Fig.9)

At the compressed air outlet (1) of the compressor, connect the pressure hose with a nut (2) and fasten clip (3). Connect the hose to the dental appliance.


Fig.9

9.3. Electrical connection

Plug the electrical cord into the mains.

The appliance is equipped with a grounded plug. Make sure this connection complies with local electrical codes. The mains voltage and frequency must comply with the data stated on the appliance label.




- Keep the socket easily accessible to ensure that in an emergency the appliance can be safely disconnected from the mains.
- Connection to the power distribution box must be max.16 A.

Electrical cable may not contact the hot parts of a compressor. Insulation could be damaged!
If any electrical cord or air hose is damaged it must be replaced immediately.

9.4. Chance in the door opening

- Disassemble the door, rectification screw (3) and the holder (2) of the hinge D (4).
- Mount holder of the hinge D to the left side of the casing.
- Rotate the door by 180°.
- Insert a spacer (5) between the hinge H(1) and the bottom side of the door.
- Mount the door.
- Disassemble the lock (6) on the door, rotate it by 180°.
- Disassemble the latch (7) on the lock, rotate it by 180°.
- Mount the lock to the door.


10. WIRING DIAGRAMS

DK50 Z, DK50-10Z, DK50 S, DK50-10S

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.


B TYPE

DK50 Z, DK50-10Z, DK50 S, DK50-10S

1/N/PE ~ 110 V 60 Hz
ELECTRIC OBJECT OF 1st CAT.

B TYPE


DK50 PLUS

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.

B TYPE


B TYPE
ST

DK50-10Z/M, DK50-10S/M

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.


B TYPE

DK50-10Z/M, DK50-10S/M

1/N/PE ~ 110 V 60 Hz
ELECTRIC OBJECT OF 1st CAT.


B TYPE

DK50 PLUS/M1a

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.


B TYPE

DK50 Plus S

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.

B TYPE


DK 50 Plus

17/N/PE-115V 60HZ
ELECTRIC OBJECT OF 1ST. CAT
TYPE B

TYPE B

ST1


DK 50 Plus/M1a

1/1/PE = 115V 60HZ
ELEKTRICKÝ PREDMET TR. I
TYP B

TYPB

ST

M	Motor of compressor
EV1	Fan of compressor
EV2	Fan of dryer
YV2	Solenoid valve of dryer
YV1	Solenoid valve of compressor
FU	Fuses 230/50-60 (T10A) 110/50-60 (T16A)
ST	Thermo switch
Cb1.Cb2	Capacitor

SP	Pressure switch
SRH	Humidity sensor
X1,X2	Terminal
X4	Connector
FA	Breaker
HL1	Glowlamp
Q	Switch

11. FIRST OPERATION

(Fig.10)

- Make sure that all stabilizers used during transport were removed.
- Check that all pressurized air line connections are secure.
- Connect to the mains.
- Start compressor at pressure switch (2) by turning switch (3) to position "I."
- For kompressor DK50 PLUS S turn the switch (29) (Fig.6) at the front part of the soundproof box to the position "I" – green light indicates that the appliance is on.

Compressor - At first operation the air tank is pressurized until it reaches a preset level when the compressor automatically switches off. As the air is used, the compressor works in automatic mode, switched on or off by the pressure switch.

Compressor with dryer - In addition, the adsorption dryer takes out the humidity from the previous compressed air during operation and it blows out the entrapped condensate via condensate outlet, which can be heard as a short hiss when compressor is stopped.

Compressor with condensation and filtration unit - Model KJF-1 filters and dehumidifies the air and automatically releases condensed liquid through the filter's discharge valve.

The compressor is not equipped with an emergency power supply.

OPERATION

In case of emergency, disconnect the compressor from the mains (pull out the mains plug).

**The compressor has hot surfaces.
Burns or fire may result if contact is made.**

During prolonged operation of the compressor, the temperature in the box may increase to over 40°C. At this point the cooling fan automatically switches on. After cooling the space to under 32°C, the ventilator switches off.

Automatic start: when pressure in the tank drops to the pressure switch's lower limit level, the compressor automatically switches on. The compressor automatically switches off after reaching the pressure switch's upper limit level.

Compressor with dryer

A correct function of the drier depends on the compressor's operation and no attendance is required. The pressure vessel need not be sludged, because the pressure air entering the air chamber is already dried.

- For correct operation of a drier it is necessary to:
- Observe the operation of the compressor in intermittent mode up to 60%. However, the time of continuous operation of the compressor should not be longer than 10 minutes.
- It is forbidden to alter the working pressures of pressure switch set by manufacturer. The operation of the compressor at working pressure lower than the switching pressure demonstrates the overload of the compressor (high air consumption) by the appliance, leakages in pneumatic distributions, failure of aggregate or drier.
- Leave the compressor connected to electric mains (do not switch off pressure switch and do not disconnect power cord) – if the humidity sensor assesses air as "insufficiently dry", it shall regenerate the drying chamber refill by intake of air from air chamber in which pressure drops and it may come to switching on the compressor several times. After the drier is regenerated, the operation of the compressor shall be automatically halted.
- Should regeneration takes longer than 1 hour and there is no air consumption, it shall be necessary to verify whether:
 - Regeneration takes place,
 - Air goes out from the outlet of solenoid valve of the drier via condensate drain,
 - There is no failure of the compressor or drier,

12. SWITCHING THE COMPRESSOR ON

(Fig.10)

Switch on the compressor at the pressure switch (2) by turning the knob (3) to position "I." (for DK50 PLUS S switch on also switch (29) Fig.6, on the front part of the compressor box), The compressor sends pressurized air to the air tank. As the compressed air is used, the pressure in the air nozzle drops to a preset level, the compressor switches on and the air nozzle fills with compressed air. After reaching the cutoff pressure the compressor turns off automatically and the cycle is repeated. Check the value of switching-on and switching-off pressure on pressure gauge. The values may be within a tolerance of $\pm 10\%$. Air pressure in air chamber must not exceed maximal permitted operation pressure.


Fig.10

Never tamper with the pressure switch (2). Adjustments are not allowed. The pressure switch (2) has been set by the manufacturer and further setting of switching on and off pressure may be carried out only by a qualified expert trained by the manufacturer.

After putting the compressor with a drier M1a to operation (after installation or several-day break) the regeneration of the drier may be started. The compressor is put into operation automatically also without pressure air take-off by the appliance, and thus increase pressure in the air jet into cut-off pressure when it turns off. Then the dryer regeneration takes place (passing of the air from air jet through dryer chamber). The air jet pressure decreases to closing pressure, the compressor turns on and increases pressure into air jet to cut-off pressure and turns off. The process of turning on and off of the compressor repeats in this way until the dryer is not regenerated sufficiently. This level is controlled in the dryer through the built-in humidity monitor – hygrostat. The regeneration process can take several minutes (5 – 15 min) – for the new or previously already regenerated dryer during the previous operation, or several tens of minutes (30 – 120 min) – for the dryer being previously in its operation „overloaded“ with water vapour (for example due to the activity regime of the compressor outside the permitted extent, during the work in the environment with high relative humidity, etc.). After finishing regeneration the whole process is stopped automatically.


Note!

The mentioned values of the regeneration times are valid for compressor activity regime at closed outlet valve, which means without pressure air take-off by the appliance. In case of air take-off the mentioned times are prolonged.

In the case the compressor does not terminate the automatic regeneration of a drier within 120 min, it shall be necessary to turn to your supplier or servicing centre.

The table states the time of regeneration cycle of the compressor with a drier M1a (under the condition that no air is taken from the compressor by an appliance).

Used jet :
DK50 PLUS 0,5mm

Compressor	Jet size	Closing pressure - cut-off pressure	Compressor activity time	Compressor break time - dryer regeneration
DK50 PLUS/M1a, DK50 PLUS S/M1a	0,5 mm	4,5 – 6,5 bar	About 60 – 70 s	about 180 – 210 s

In case of deviation from described activity please verify:

Connection of the solenoid valve – compare arrow orientation on the valve body, which also marks correct air flow direction when regenerating

Possible leakage via compressor untightness – close the outlet valve on the compressor, turn on the compressor and leave it in operation until it turns off at cut-off pressure. Turn off the circuit breaker and watch pressure on manometer, which means the pressure in air jet. The pressure decrease can not be bigger than 0.2 bar in 2 hours.

The use of the correct jet – a jet is mounted between solenoid valve and opening (see Table). When assembling it is important to keep jet orientation – aim the surface with a screwed cone at the valve.

MAINTENANCE

13. MAINTENANCE SCHEDULE

Maintenance that must be performed	Chapter	Time interval	Performed by
- Release condensate - Compressor without air drier - At high air humidity - Compressors with air drier - Compressors with condensation unit : - from filter - from pressure vessel	14.1	1 x week 1 x day 1 x week, check function 1 x week, check function 1 x week	operating staff operating staff operating staff operating staff operating staff
Check safety valve	14.2	1 x year	qualified technician
Replacement of input filter	14.3	1 x 4 years or after 8000 hours	qualified technician
Replacement of output filter in dryer	14.4	1 x 2 years	qualified technician
Replacement of filter in condensation unit	14.5	1 x year	qualified technician
- Check tightness of joints - Overall examination of device	Service documentation	1 x year	qualified technician

14. MAINTENANCE

Repair work beyond normal maintenance can be performed only by qualified personnel or the manufacturer's representative.
Use only spareparts and accessories approved by the manufacturer.

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).


TO ENSURE THAT THE COMPRESSOR WORKS CORRECTLY, PERFORM THE FOLLOWING MAINTENANCE TASKS AT REGULAR INTERVALS (CHAPTER 13):

14.1. Condensation drain valve

Compressors (Fig.11, Fig.12)

During regular use, release condensation from the pressure tank. Switch off the compressor at the mains. Reduce air pressure in the appliance to max. 1 bar by releasing air via a connected device. Place the hose with the drain valve into a container prepared in advance (for compressors DK50 PLUS place the vessel under release valve) and open the drain valve (1). Wait until condensation is fully drained from the pressure tank. Close drain valve (1).

Fig.11

DK50

Compressors with condensation and filtration unit (Fig.15)

During regular use, condensation is automatically released via the release valve of the condensation unit filter. To check that the automatic drain is working properly, open the valve (4) of the drain vessel (2) by turning to the left. Release a small amount of condensate from the vessel. Close the valve (4) by turning to the right.

Compressors with air dryer

In the case of a regular operation condensate is automatically excreted via air dryer and it is entrapped in a bottle. Take out the bottle from a holder, release blowing-out suppressor and pour out the condensate.

Fig.12

DK50 PLUS

If necessary, it is possible to connect the set for condensate discharge onto the condensate outlet (see Chap. PARTS LIST - Auxiliary Equipment).

For compressor models DK50 S, DK50-10 S and DK50-10S/M the case must be removed before beginning the following procedures.

For DK5010-S/M - first disconnect the regeneration valve (Fig.5).

For DK50 PLUS S, DK50 PLUS S/M - unlock the lock on the door and open the housing door (Fig.6)

14.2. Safety valve check

(Fig.10)

When the compressor is operated for the first time, make sure that the safety valve is working properly. Turn screw (4) of safety valve (1) several rotations to the left until the safety valve releases air. Let the safety valve blow out for only a few seconds. Turn screw (4) to the right until it seats, closing the valve.

The safety valve must never be used for depressurizing the air tank. It could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer. Adjustments are not permitted.

Warning! Compressed air can be dangerous. Wear eye protection when blowing air out.

14.3. Replacement of input filter

(Fig.13)

It is necessary to replace the input filter (1) located in the lid of crank box of a compressor.

- Pull out the rubber plug (2) using a hand.
- Take out the used and contaminated filter.
- Insert new filter and put on a rubber plug.


Fig.13

14.4. Replacement of output filter in air dryer

Before beginning, reduce the air pressure in the tank to zero and disconnect the appliance from the mains.

(Fig.14)

In the case of a regular operation of a dryer it is necessary to replace the dryer filter in the upper part.

- Remove the dryer plug (1) by unscrewing to the left from the dryer head (4).
- Put in a new filter (2).
- Refit the dryer plug to the dryer head (4) and tighten it to the right.


Fig.14

14.5. Replacement of filter in condensation and filtration unit

Before beginning, depressurize the air tank to zero and disconnect the appliance from the mains.

(Fig.15)

In the case of a regular operation of a condensation unit it is necessary to replace the filter inside the filter with automatic desludging.

- Release a safety lock (1) on the filter vessel by its pulling downwards, slightly rotate the filter cover (2) to the left and take it out.
- Unscrew the filter holder (3) by its rotation to the left.
- Replace the filter and fix the new one by rotation of the holder to the right back on the filter body.
- Replace the filter cover and secure it by turning to the right until the safety pin locks.


Fig.15

15. STORAGE

If the compressor will not be used for a prolonged time period, drain any condensate from the air tank. Then turn on the compressor for 10 minutes, keeping the drain valve open (1) (Fig.11, Fig.12). Switch off the compressor by switch (3) at pressure switch (2) (Fig.10), close the drain valve and disconnect the appliance from the mains.

16. DISPOSING OF THE APPLIANCE

- Disconnect the appliance from the mains.
- Release air pressure in the pressure tank by opening the drain valve (1) (Fig.11, Fig.12).
- The components of the product are non-toxic.
- Dispose of the appliance following all environmental regulations.

17. REPAIR SERVICE

Guaranteed and post-guarantee repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

The manufacturer reserves the right to make changes to the appliance without notice. Any changes made will not affect the functional properties of the appliance.

18. SOLVING PROBLEMS

Caution! Before proceeding, depressurize the air tank to zero and disconnect the appliance from the mains.

After troubleshooting and reassembly of the drier, it is necessary to carry out the drier regeneration by starting the compressor with a small air offtake (compressor working mode 20-30%) up to the automatic stop of regeneration.

Troubleshooting can be performed only by qualified personnel.

FAILURE	POSSIBLE CAUSE	REMEDY
Compressor does not start	No voltage in pressure switch Disconnected winding of motor, damaged thermal protection Faulty capacitor Seizure of piston or another rotary part Pressure switch does not switch on	Check voltage in socket Check fuse – replace faulty one Loosen terminal – tighten it Check power cord – replace faulty one Replace motor or re-wind it Replace capacitor Replace damaged parts Check the function of pressure switch
Compressor often switches on	Air leak in pneumatic distribution system Leaking check valve Greater volume of condensed liquid in pressure vessel	Check pneumatic distribution system – seal loose joint Clean valve, replace seals, replace valve Drain condensed liquid
Prolonged running of compressor	Air leak in pneumatic distribution system Worn piston ring Contaminated input filter Dirty filter in the dryer Defective solenoid valve	Check pneumatic distribution system – seal loose joint Replace worn piston ring Replace contaminated filter with the new one Change outlet filter and inspect dessicant Repair or change the valve
Compressor is noisy (knocking, metal noises)	Damaged bearing of piston, piston rod, motor bearing Loose or cracked spring	Replace damaged bearing Replace damaged spring
Dryer doesn't dry (condensed water in the tank)	low operating pressure	reduce the air offtake, check the source efficiency, check possible leaks in the distribution
	inoperative solenoid valve	repair or change the valve
	clogged nozzle of the regenerating air	clean or replace the nozzle, use correct size of the nozzle (see Product maintenance)
	inoperative cooler ventilator	replace ventilator check supply of electric energy
	leak of white liquid through the solenoid valve	dismantle the chamber, replace drying substance, bottom filter and sealing and inspect tightness, apply soapy water onto O-rings of the nuts
Drying is noisy	low operating pressure	reduce the air offtake, check the source efficiency, check possible leaks in the distribution
	defective solenoid valve	replace the valve
	damaged damping substance in the condensate vessel	replace damping substance or the vessel
	damaged pressure hose	replace pressure hose

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CERTIFICATE:

This product is made according to strict manufacturing standards and meets all technical specifications.

GUARANTEE:

This product is guaranteed by the manufacturer for 24 months from the date of sale providing that all terms of the guarantee are met.

- To request repairs under this guarantee, submit a claim form along with the product's quality certificate together with the product itself.
- During the guarantee period, any defect caused by the manufacturer will be corrected at no charge.
- This guarantee does not cover damage during transport or damage caused by incorrect storage.
- This guarantee does not apply to damage resulting from incorrect operation.
- This guarantee does not apply if the product is used for any purpose other than those described in the user's manual or explicitly approved by the manufacturer.
- This guarantee does not apply if the product has been adjusted or altered contrary to the manufacturer's instructions.
- This guarantee does not apply if the product was obviously incomplete or defective at the time of purchase and the buyer knowingly accepted it.

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